

Using Your IN-APP REPORTS





To access your In-App Reports, go to Main Menu > Reports


Your Eversense 365 CGM System automatically creates 4 reports you can view at any time — right in the app. These reports can help you and your healthcare provider (HCP) fine tune your diabetes management plan.

Glucose Pie Chart report:

- Shows time spent within, below, and above your established target range.
- Options for either 1, 7, 14, 30, or 90 days of data display.
- Comparing results from these different intervals may help you and your HCP see how your diabetes care plan is working, and perhaps make adjustments.

Time in Range report:

- Displays percentage of time you've spent in pre-set glucose ranges, based on standardized guidelines.
- Tapping on the information icon displays more information about these fixed ranges.
- Options for either 1, 7, 14, 30, or 90 days of data.



Weekly Modal Summary Report:

- Provides a seven-day summary of your glucose information in a 24-hour line graph.
- Easy-to-see glucose readings at different times of day.
- Red bars show range of glucose variability at a given time.
- Scrolling down shows additional glucose statistics at pre-set times.

Glucose Statistics Report:

- Shows average, low, and high glucose values based on preset times.
- Options for either 1, 7, 14, 30, or 90 days of data.
- Comparing results from these different intervals may help you and your HCP see how your diabetes care plan is working, and perhaps make adjustments.



All reports can be emailed directly through the app.

© 2024 Senseonics, Inc. All rights reserved. © 2024 Ascensia Diabetes Care Holdings AG. All rights reserved.



Distributed by Ascensia Diabetes Care 5 Wood Hollow Road Parsippany, NJ 07054 844-SENSE4U (736-7348) eversensecgm.com

Manufactured by Senseonics, Inc. 20451 Seneca Meadows Parkway Germantown, MD 20876-7005 USA The Eversense' 365 Continuous Glucose Monitoring (CGM) System is indicated for continually measuring glucose levels for up to one year in people (18 years and older) with diabetes. The system is indicated for use to replace fingerstick blood glucose (BG) measurements for diabetes treatment decisions. Fingerstick BG measurements are required for calibration one time a week after day 13, and when symptoms do not match CGM information or when taking medications of the tetracycline class. The sensor insertion and removal procedures are performed by a health care provider. The Eversense 365 CGM System is a prescription device; patients should talk to their health care provider to learn more.

For safety information, see bit.ly/eversensesafety

Eversense, Eversense 365 Continuous Glucose Monitoring, and the Eversense logo are trademarks of Senseonics, Incorporated. Other brands and their products are trademarks or registered trademarks of their respective holders and are used solely for informative purposes. No relationship or endorsement should be inferred or implied. PP-SENS-365-GBL-0053